

# The Importance of Factorial Design in Tissue Engineering and Biomaterials Science: Optimization of Cell Seeding Efficiency on Dermal Scaffolds as a Case Study



Alexandra Levin<sup>1</sup>, Vaibhav Sharma<sup>1</sup>, Lilian Hook<sup>2</sup>, Elena García-Gareta<sup>1</sup>

<sup>1</sup>Regenerative Biomaterials Group, RAFT Institute, Mount Vernon Hospital, Northwood, UK

<sup>2</sup>Smart Matrix Ltd, Leopold Muller Building, Mount Vernon Hospital, Northwood, UK

**Background:** We present a case study to show the usefulness and importance of using experimental factorial designs in tissue engineering and biomaterials science.

We used a full factorial design (2x2x2x3) to solve a routine query in biomaterials research: the optimisation of cell seeding efficiency for pre-clinical *in vitro* cell studies, the importance of which is often overlooked. Moreover, tissue-engineered scaffolds can be cellularised to form implantable tissue constructs, where the cell seeding method must be reliable and robust.

**Aim:** to optimise cell seeding efficiency on dermal scaffolds for *in vitro* pre-clinical studies using full factorial design

## Method & Results

**Dermal scaffolds:** Integra® (collagen/GAG+silicone) and bovine Smart Matrix® (fibrin/alginate) cut to 6mm diameter.

- Homogenous structures of open, interconnected macro and micro-pores (A,B).
- Nano-pores and densely packed nano-fibres only observed for Smart Matrix® (B).
- Integra® is mechanically stronger and retains larger volumes of liquid (~125 µL) than Smart Matrix® (~25 µL) due to the silicone backing layer (C).


**Cells:** primary normal human dermal fibroblasts (pnhDFs) from a single donor.

- Cells maintained their spindle-shaped morphology throughout the study (A).
- Immunostaining of cells for Ki67 (red) and actin (green) suggests cells were proliferative at the time of the experiments (A).
- Cells were viable at the time of the experiments (B).



**Experimental design:** Variables and levels investigated (see matrix below), based on our experience with these materials, were:

- 1) cell passage number (5 or 10)
- 2) cell seeding density (1.25x10<sup>5</sup>, 2.5x10<sup>5</sup> or 5x10<sup>5</sup> cells in 200 µL)
- 3) scaffold disc to well plate surface area ratio (1:1 or 1:6)
- 4) attachment incubation time (3 h or 24 h).

- Full factorial experimental design (2x2x2x3).
- For each individual set of experimental conditions n=3.
- Cell seeding efficiency was quantitatively assessed using alamarBlue®, a metabolic redox assay and calculated as percentage of cells remaining on the scaffolds.
- A standard curve was created for each passage number and attachment incubation time.
- Cell seeding was qualitatively assessed by histological processing and microscopy.



| Attachment | Integra®       |  |                |  | Cell Seeding Density | Smart Matrix®  |        |                |  | Incubation |          |    |
|------------|----------------|--|----------------|--|----------------------|----------------|--------|----------------|--|------------|----------|----|
|            | Passage Number |  | Passage Number |  |                      | Passage Number |        | Passage Number |  |            |          |    |
| 3          | P5             |  | P10            |  | 125000               | P5             |        | P10            |  | 3          |          |    |
| 24         |                |  |                |  |                      | 250000         |        |                |  |            | 24       |    |
| 3          |                |  |                |  |                      |                | 500000 |                |  |            |          | 3  |
| 24         |                |  |                |  |                      |                |        |                |  |            |          | 24 |
| 3          |                |  |                |  |                      |                |        |                |  | 3          |          |    |
| 24         |                |  |                |  |                      |                |        |                |  | 24         |          |    |
| Time (h)   |                |  |                |  |                      |                |        |                |  |            | Time (h) |    |



**Main effects and interactions:** attachment incubation time had a strong negative effect on seeding efficiency (B).






- Scaffold disc to well plate surface area ratio had a strong negative effect for Smart Matrix.
- Increasing passage number and cell seeding density had a negative effect for both scaffolds
- Statistical analysis showed that for Smart Matrix® the main effects of attachment incubation time and scaffold disc to well plate surface area ratio were statistically significant and so was their interaction (C).
- For Integra® only the main effect of attachment incubation time was statistically significant.

**Microscopy:** Phase-contrast light microscopy of empty wells revealed a ring of cells left behind following scaffold removal from 24 well plates (A).

- Fewer cells were left behind in 96 well plates: using a 96 well plate restricts cell seeding adhesion to the scaffold.
- H&E staining of seeded scaffolds revealed a layer of cells at the top of the scaffold where they were seeded (B).
- Qualitatively fewer cells were observed as the seeding efficiency decreased.

**Visual representation of results:** in order to more clearly observe the effects and interactions of the different variables and find the optimum combinations that should be used for each scaffold, we propose 2 different visual representations of the data:

- 1) 3D graphs (A),
  - 2) the matrix above, was filled with results and a colour key was assigned to values (B).
- For Integra®, highest efficiencies were found when:
- 1) 5x10<sup>5</sup> cells at P5 where seeded on scaffolds placed in 96 well plates (1:1) and incubated for 3h (60.2%),
  - 2) 1.25 x10<sup>5</sup> cells at P10 were seeded on scaffolds placed in 24 well plates (1:6) and incubated for 3h (59.1%).
- For Smart Matrix®, highest efficiencies were found when:
- 1) 1.25x10<sup>5</sup> cells at P10 where seeded on scaffolds placed in 96 well plates (1:1) and incubated for 3h (105.3%),
  - 2) 1.25 x10<sup>5</sup> cells at P5 were seeded on scaffolds placed in 96 well plates (1:1) and incubated for 3h (91.1%).



## Conclusion:

- Our study design could save time and resources and the optimum seeding conditions should be investigated for individual scaffolds.
- Our study can be easily translated to other cell types and biomaterials, where multiple interacting variables can be thoroughly investigated for better understanding cell-biomaterial interactions.

## References & Acknowledgements:

Levin *et al.* 2018. J Tissue Eng 9:1-14; Sharma *et al.* 2016. Biomed Mater 11:055001. Work supported by the Restoration of Appearance and Function Trust (UK, charity No. 299811).

